

MAINTAINING YOUR EQUIPMENT

- **MASK INTERFACE:**
 - Daily: Wipe down cushion/pillow with damp washcloth or CPAP wipes
 - Weekly: Disassemble mask and wash cushion/pillow with gentle soap and warm water
 - Monthly: Replace cushion/pillow
 - Replace entire mask including frame and headgear every 3 months

- **TUBING:**
 - Weekly: Detach from machine and rinse with soap and water; hang to dry
 - Replace tubing every 3 months

- **WATER CHAMBER:**
 - Daily: Empty and refill with distilled water
 - Weekly: Wash thoroughly with soap. Be sure to rinse completely and let dry

- **NON-DISPOSABLE FILTER (DARK BLUE):**
 - Do not wash
 - Replace non-disposable filter every 3 months

- **DISPOSABLE FILTER (LIGHT BLUE):**
 - Do not wash
 - Replace every 2-4 weeks as discoloration appears